# GARGOLA HEALTHCARE Hospital supplies and services

# Imported by Gárgola Healthcare, S.L.

# National Number 205428.8 **General Council of Pharmaceutical Colleges**





**AUTO TEST BOSON** 

Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

REF 1N40C5-2

LOT 21051707

2022-11



Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian, 361021, P.R.China. www.bosonbio.com



Importer: Gárgola Healthcare, S.L. Calle Vázquez de Mella 6, 2A, 50.009 Tax number: B99498701



Rapid SARS-CoV-2 Antigen Test Card
















CN 205428.8

COVID-19





NO

iagnostisk selvtest for hurtig påvisning av SARS-CoV-2





Advarse



(1) Uwaga



Atenção

1 Attention

ES



SSV









# TESTSEALABS

### Covid-19 Rapid Test Kit

Chek

### **Testsealabs Covid Antigen Test Cassette**

### Designed for the qualitative detection of SARS-CoV-2 antigen on swabs

The Covid-19 Antigen Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of COVID-19 antigen in a nasopharyngeal swab sample to aid in the diagnosis of COVID-19 viral infection.

### **Product main features:**

- Test time: in 15 minutes / Sensitivity: 92.1%; 95% CI: 87.7% 93.9%) / Specificity: 98.3%; 95% CI: (8;94.2% - 100%) / Sample: swab / Storage state: 4 - 30°C
- Shelf life: 2 years / Available formats: Box of an individual Test, Box of 20 Tests.

One individual test box content: 1 test cassette, 1 sterile swab, 1 extraction tube, 1 product leaflet.

Contained in a box of 20 units: 20 test cassettes, 20 sterile swabs, 20 extraction tubes, 20 product inserts, 20 extraction buffers.



GICA

## Testseaelabs ® (COVID-19SARS-CoV-2) Antigen Schnelltest (BfArM GZ: 5640-S-214/21)

### **Example type:**

Anterior nasal swab, non-invasive and painless sampling;

- Easy to operate, suitable for COVID-19 family and individual self-test;
- Self-test, continuous self-monitoring of health status, early detection, early isolation, early treatment
- High sensitivity, forms the first line of defense against COVID-19 detection, reduces pressure from medical institutions

Art. Ref: PBLCATK05

- Suitable for multiple scenarios:

Tests in medical institutions; screening before resuming work and school, continuous monitoring, etc.

- All necessary reagent is provided and no equipment is required;
- Time-saving procedures, results are available in 15 minutes;
- Storage temperature: 4 ~ 30 °C No need for cold chain transport;
- Specification: 20 tests / box, 5 tests / box, 1 test / box.





## **Professional use test** Home use test



